

DaVinci Feed Tube Set

DaVinci Feed Tube Set - 70 Series - 0.5 ml Catlog No.: 70-02-23-29




Figure:1

Each 0.5 ml DaVinci feed tube set syringe kit includes Metal Body, Luer lock nozzle, Plunger, Glass barrel, Integrated valve system, Feed Tube, Sinker valve, Feed cannula, Spare kit (Springs, Valve balls & Washers), Detailed operating manual.

Specifications:

1. Plunger tightnes without 'O' Ring.
2. Barrel made of Borosilicate glass
3. Clear visibility of Dose
4. Barrel ID 5.6 mm \pm 0.1mm
Length 46.25 mm + 0.10 mm
5. Plunger made of Stainless Steel.
6. HS Code: 90183100

Table 1

Parts in Contact with Liquid For FTS			
S.No.	Component	Part Number	Material
1	Feed Tube	23	Silicone
2	Glass Barrel	7	Borosilicate glass
3	Washer, Barrel	15	FPM
4	Plunger	6	Stainless Steel
5	Valve System	11	Ni-Cr Brass
6	Washer, Valve	13	PTFE
7	Spring	12	Stainless steel
8	Ball Valve	17	Stainless steel
9	Sinker	24	Ni-Cr Brass
10	Feed cannula	26	Ni-Cr Brass




Figure:2

TOP SYRINGE
MANUFACTURING COMPANY PRIVATE LIMITED

10/11, Prospect Chambers Annexe, 317-21, Dr. D.N. Road, Fort, Mumbai - 400 001, Maharashtra, India.

Tel.: +91-22-6520 8849/50, +91-22-6520 6102/03 Fax: +91-22-2204 5352 • www.glass-syringe.com • sales@glass-syringe.com